Grazedean Ltd

ARC PLUS - ARGON UNIT for BOWA HF GENERATORS



AND ARC CART



Completing your work station



Tel +44 (0)1202 584645 Fax +44 (0)1202 593824

email sales@grazedean.com www.grazedean.com

ARC PLUS - ARGON UNIT





With the ARC PLUS argon module and accessories arcs can be generated at **less than 10 watts**, reducing the risks in fine tissue areas.

| Bronchoscopy | Argon Flex |
|------------------|--|
| Gastroenterology | Argon Flex, GastroCUT (Polypectomy, Papillotomy) |
| General surgery | Argon Ligation, Argon Open, Ligation, Standard |
| Liver Surgery | Argon open, Argon ligation |

Disposable argon flexible probes are available as well as the standard reusable. Enhanced design ensures consistent ARC generation at low wattages is achieved. More information is available - for a full catalogue please contact us.



Argon flexible probes are available for BOWA and other manufacturers machines in lengths from 1.5m to 3m long and 1.5 - 3.2mm Ø available. The rounded ceramic tip gives excellent feeding properties into a working channel and long life. Argon rigid hand control probes are availble with needles or blades in various lengths, Long life is achieved using with ceramic inserts.



Grazedean Ltd

Tel +44 (0)1202 584645 Fax +44 (0)1202 593824 email sales@grazedean.com www.grazedean.com

Technical specification for ARC PLUS

BOWA ARC PLUS can be used in all areas of HF surgery in combination with all ARC generators. It can be located both on ceiling supply units and also on a mobile ARC CART. The communications interface provided continuously monitors the interaction with the ARC generator. The unit complies with the currently required applicable standards and can be maintained by the hospital technician himself without special measuring equipment. Fault memories can be read out without difficulty. The technician can update software simply by changing a board and subsequently calibrating the unit.

| Gas supply | |
|-------------------------------------|---------------------------|
| Gas type Argon Purity 4.8 or better | Class 2.1 A |
| Adjustment range 0,1 to 1,0 l/min | Adjustment step 0,1 I/min |
| Adjustment range 1,0 to 3,0 l/min | Adjustment step 0,2 l/min |
| Adjustment range 3,0 to 9,5 l/min | Adjustment step 0,5 l/min |
| Gas input pressure between | 2,5 and 4,5 bar |
| Max. gas output pressure | 2 bar |

Two gas tank connections, automatic or manual changeover of gas tanks and optical level indicator

PURGE rinsing function with electrode dependent gas pressure control

Safety devices

Continuous flow monitoring, fault indication with display of error number, hose blockage detection, continuous self-test, warning tone, continuous pressure monitoring at distal end

Documentation

Recording of fault conditions and operator errors by storing within the unit Data recall facility

External outputs Optical interface for communication with ARC HF Surgery Generators

| Mains connection | |
|---------------------------------|-------------------|
| Power consumption | 40 VA |
| Mains frequency | 50 / 60 Hz |
| Equipotential bonding connector | Yes |
| Voltage range | 85 bis 264 V A.C. |

| Dimensions and weights | | |
|--|--|--|
| External dimensions (W x H x D) | 430 x 94 x 400 mm | |
| Net weight ca. | 7,3 kg | |
| Protection class to EN 60601-1 | 1 | |
| Type according to EN 60601-1 | CF in conjunction with ARC HF Generators | |
| Classification according to EC Directive 93/42/EEC | Class IIa | |




Tel +44 (0)1202 584645 Fax +44 (0)1202 593824 email sales@grazedean.com www.grazedean.com

ARC CART

The ARC HF generators and ARC PLUS are often combined to form an unbeatable all-in-one package, which can be conveniently located on the BOWA ARC CART or on ceiling supply units.

The ARC CART can be configured individually, e.g. shelves, drawer, basket, gas tank holder.

The whole range of foot switches has ergonomically arranged pedals.



Also available tissue sealing devices and smoke evacuation.

